

## Instructions for Checking and Tightening Screws:

1. **Inspect Regularly:** Set a routine to inspect the screws on your optical frame at least **once a month**. Pay close attention to the hinge areas and nose pads, as these are common areas for screws to loosen.
2. **Use Proper Tools:** Use the Uni tool provided in the case. Ensure that the components of the device are fitted and adjusted correctly before use. Tighten the hinge screws as appropriate to maintain the stability of the frame.
3. **Gently Tighten:** Carefully tighten any loose screws by turning them clockwise until snug. Be cautious not to overtighten, as this can strip the screw or damage the frame.
4. **Seek Professional Help if Needed:** If you're unsure about tightening screws or if you encounter any issues during inspection, visit an optician or eyewear professional for assistance.

By following these simple steps and incorporating regular maintenance into your eyewear care routine, you can ensure the longevity, comfort, and performance of your optical frames.

